

MAICO Audiometer MA 25

Light and independent for basic screening

Sleek and Small Design

The super-lightweight MA 25 weighs in at only 980 g. Small enough to hold in the palm of your hand! Tucked into its carrying case, the MA 25 is one of our smallest, lightest portables.

Accurate and Dependable Engineering

Because it is portable, the MA 25 has been engineered to be tough and reliable. Take the MA 25 everywhere with confidence. User-friendly knobs make adjusting frequency and hearing levels easy and comfortable!

Full Range of Test Options

The MA 25 offers a full range of testing opportunities. Frequencies range from 125 Hz to 8000 Hz and signal intensities run from $-10 \text{ dB}_{\text{HL}}$ to $100 \text{ dB}_{\text{HL}}$. Sinus and warble tone (pulsed and continuous) make it complete.

Features at a glance

- Air Conduction with 11 frequencies from 125 to 8000 Hz
- Customize the frequencies to your own needs
- Can be operated with 3 standard AA batteries or eco-friendly with a medically approved external power supply (optional)
- Pure, pulse, warble tones
- Exceptionally lightweight, only 980 g
- MA 25e extends functionality of MA 25



Technical Data

MA 25 SPECIFICATIONS

Test Signals	Sinus and warble tone (pulsed and continuous)
Test Frequencies Hz	125, 250, 500, 750, 1000, 1500, 2000, 3000, 4000, 6000, 8000
Level Steps	5 or 1 dB
Hearing Level Range with Headset DD45	-10 dB _{HL} - 100 dB _{HL} 70 dB _{HL} at 0.125 kHz, 90 dB _{HL} at 0.25 and 8 kHz
Audiometer Class	Class 4 according to EN 60 645-1

GENERAL

Dimensions	22.5 cm x 18 cm x 5.5 cm
Weight	980 g
Power Supply	3 AA batteries (100 hours life time typically) or optional power supply

STANDARD COMPONENTS

- Audiogram Pad
- Operators Manual
- DD45 Headset Assembly

OPTIONAL ACCESSORIES

- Hughson-Westlake test
- Patient response switch
- Power supply
- Carrying case
- Headset Holmco 8103, DD65 and TDH 39

STANDARDS

EN 60 601-1
EN 60 645-1, class 4,
according to medical device directive 93/42/EEC



MA 25E EXTENDS FUNCTIONALITY OF MA 25

The MA 25e extends the MA 25 functionalities with an automatic test according to Hughson-Westlake. Parameters can easily be adapted to your personal preferences. The Talk Forward function makes the MA 25e easy to work with, particularly in sound booth installations.

Standard Components				
	MA 25 device	Headset DD45	Audiogram pad	
	Optional Accessories			
		Headset Holmco 8103	Headset DD65	Headset TDH39
				
Patient response switch				

MAICO Diagnostics GmbH

Sickingenstr. 70-71 · 10553 Berlin · Germany
Tel.: +49 30/70 71 46-50 · Fax: +49 30/70 71 46-99
sales@maico.biz · www.maico.biz