

David Smith, US Olympic Team volleyball player and Oticon Dynamo wearer, enjoys a day in the park with his wife, Kelli, and their son, Cohen.



Visit oticonusa.com to hear David and Kelli's story about Dynamo.

Oticon | Dynamo

People First

People First is our promise to empower people to communicate freely, interact naturally and participate actively

The power of understanding



10592 15555-0070/10.15

“ It's been great just to be able to perceive sounds that I've never perceived before because they were outside of my range. And the power that Dynamo gives me just makes the sound so much easier to understand.

”
– David Smith



The freedom to hear conversation like never before.


www.oticonusa.com

oticon
PEOPLE FIRST



oticon
PEOPLE FIRST

Imagine being able to hear more of the details of speech than you have in the past, with a super power hearing device that's programmed for you and only you. And to be able to understand more of the nuances of everyday conversations like never before.

A powerful advance in hearing technology.




Innovations like this don't come along every day. That's why you owe it to yourself to take a closer look at Dynamo by Oticon. Our newest and most innovative super power hearing instrument to-date, Dynamo not only features our most advanced technology, it gives you our most power ever in a compact BTE13 Super Power device.



Connect with a world of entertainment, information and ideas.

Now your phone calls, favorite music, mobile and computer video chats, podcasts, TV, live performances and presentations can be streamed into your hearing devices via the Streamer Pro. And it can be controlled with your smartphone with our free ConnectLine App. ConnectLine is designed to give you an easier way to connect with the world of entertainment, information and ideas.



Introducing Dynamo by Oticon

Using your brainpower to understand more.

Hearing actually occurs in your brain, not your ears. And because your brain is different from everyone else's, you hear differently too. Dynamo harnesses the power of Oticon's innovative BrainHearing™ technology. This allows Dynamo to be customized to your unique hearing profile and personal sound preferences. Then it delivers sound in the way your brain is best able to understand it.

BrainHearing™ technology:

Helps both ears work together

Recognizes and preserves natural speech

Protects the clarity of sound

Coordinates how sound is best delivered to your brain

Speech Guard E keeps you in the conversation.

Different listening environments pose different challenges, but Dynamo levels the playing field. Speech Guard E works in partnership with Speech Rescue to protect the clarity of sound. In stable sound environments, the system responds slowly, and simply makes sound clear and undistorted, providing your brain with the best signal possible. When the sound environment changes abruptly, the volume adjusts swiftly and automatically, to keep the sound from becoming too loud and uncomfortable or disappearing unexpectedly. This results in a clear, stable and intact signal that the brain can decode with as little listening effort as possible.

With Speech Rescue™ you'll hear with more clarity.

Speech Rescue is a unique sound processing feature that captures spoken sound and makes it more understandable. Other power hearing instruments compress sounds in the hearing spectrum. Dynamo with Speech Rescue actually rescues the important speech frequencies that you might otherwise lose. Subtle sounds like "s" and "th" are then repositioned within your individual range of hearing to give you greater access to speech sounds, so you'll be able to better enjoy the rhythm and flow of every conversation.



Dynamo is available in seven colors that complement natural skin and hair tones.

Consult with your hearing care professional today to hear more about the difference Oticon Dynamo can make in your life.

“ It’s been great just to be able to perceive sounds that I’ve never perceived before because they were outside of my range. And the power that Dynamo gives me just makes the sound so much easier to understand.

”

– David Smith

David Smith, US Olympic Team volleyball player and Oticon Dynamo wearer, enjoys a day in the park with his wife, Kelli, and their son, Cohen.

Visit oticonusa.com to hear David and Kelli’s story about Dynamo.



The freedom to hear conversation like never before.

Imagine being able to hear more of the details of speech than you have in the past, with a super power hearing device that’s programmed for you and only you. And to be able to understand more of the nuances of everyday conversations like never before.